

# Eclipse® Accessories

Personal Ambulatory Oxygen System™ (PAOS)



## **Mobilize.**

Go everywhere and anywhere

## **Accessorize.**

Make traveling easy and convenient

## **Personalize.**

Enhance your Eclipse to complete your active lifestyle

**Travel to all of life's destinations. Whether  
around your home, your town or your own world.  
Surpass and outshine all expectations.**



 **SEQUAL®**  
The Smart Choice  
[www.sequal.com](http://www.sequal.com)

**MOBILIZE. ACCESSORIES. PERSONALIZE.**

Go to your favorite places around the globe. Travel with ease around your home, your town or the world.

The Eclipse<sup>®</sup> 2 Personal Ambulatory Oxygen System™ (PAOS) with AutoSAT™ Technology was developed with the patient in mind. The Eclipse is the only oxygen concentrator on the market that provides both continuous flow and pulse dose options in a 24/7 portable package. This allows oxygen-dependent patients the freedom of mobility and travel, even extended overnight visits. The Eclipse is truly the only, "single solution" oxygen system.

Eclipse accessories are available from your oxygen provider. Add a Protective Cover to shield rain from your Eclipse or use the Travel Case to store all of your power accessories. Our accessories were created to aid in achieving an active lifestyle deserved by all.

**TRAVEL AND MOBILITY ACCESSORIES**

**Eclipse Protective Cover Item #5052**

When zipped securely around the Eclipse, this protective cover allows limited water resistance and protection against vertically falling drops of water when tilted up to 15° from vertical (IPX2 rating)

Reduces risk of external cosmetic damage to the Eclipse unit  
Clear plastic window allows unobstructed viewing of the Eclipse control pane



**TRAVEL AND MOBILITY ACCESSORIES**

**Eclipse Travel Case Item #5010**

Cushioned interior comes with Velcro adjustable, padded, and removable partitions  
Clear zippered plastic pocket conveniently stores paperwork or small personal items  
External Velcro flap on back aids fitting over the Eclipse Profile Cart (#7092) and the Eclipse Large Wheel Cart (#7090)



## TRAVEL AND MOBILITY ACCESSORIES

### Eclipse Travel Accessory Kit Item #5093

Contents include Eclipse Travel Case (#5010), Eclipse Power Cartridge (#7082), Desktop Charger with Power Adapters (#7112), Eclipse 2 DC Power Supply with Vehicle Accessory Outlet Plug (#4142), and Eclipse Quick Start Guide



## TRAVEL AND MOBILITY ACCESSORIES

### Eclipse Wheelchair Pack Item #5220

Hangs on wheelchairs with push handles, and attaches using adjustable buckle style clips. The Eclipse fits neatly into the main compartment and is secured with compression straps and a top flap that clips shut. Clear plastic window on top flap allows unobstructed viewing of the Eclipse control panel.

## TRAVEL AND MOBILITY ACCESSORIES

### Eclipse Large Wheel Cart with Telescoping Handle Item #7090

Six inch wheels help facilitate stability and maneuverability over challenging terrain. Handle slides up and down with different settings for personalized adjustment.



## TRAVEL AND MOBILITY ACCESSORIES

### Eclipse Scooter Pack Item #5221

Soft mesh pocket easily slips over the back of scooter seat for a comfortable and sturdy fit. The Eclipse fits neatly into the main compartment and is secured with compression straps and a top flap that clips shut.

**All Travel and Mobility Accessories are interchangeable between Eclipse and Eclipse 2**

### HUMIDIFICATION

#### Eclipse Humidifier Adapter Kit Item #7116

Supports use of a humidifier with Eclipse and Eclipse 2  
Attaches humidifier bottle to Eclipse unit;  
humidifier bottle not included  
Interchangeable between Eclipse and Eclipse 2



### POWER ACCESSORIES

#### Eclipse Power Cartridge, 195 Watt-hour Item #7082

Provides 1.3 to 4.4 hours of power depending on flow rate and settings  
Recharges using AC or DC power,  
or Eclipse Desktop Charger #7112  
Interchangeable between Eclipse and Eclipse 2  
Comes standard with each Eclipse unit



### POWER ACCESSORIES

#### Desktop Charger with Power Adapters Item #7112

Charges the Eclipse Power Cartridge in less than 3 hours  
Use anywhere in the world with 100 VAC to 240 VAC, and 50-60 Hz  
Comes with detachable North American Power Cord #3890 for standard 120 VAC household outlets  
North American Power Cord #3890 can be replaced with other SeQual-authorized power cords, such as European, UK, and Hospital grade power cords



FOR ADDITIONAL ACCESSORIES, VISIT [WWW.SEQUAL.COM](http://WWW.SEQUAL.COM)

#### North America

SeQual Technologies Inc.  
(858) 202-3100  
(800) 826-4610 (US & Canada Only)

#### Europe

SeQual Technologies Europe B.V.  
+31(0) 40 290 25 60

#### Asia & Pacific Rim

SeQual Technologies Co., Ltd.  
886-2-2700-2737

