

# SIZE AND FITTING TEMPLATE




## SeQual's solution for a better fitting CPAP mask



Shown in blue for display only.

1. Print out this sheet.
2. Use scissors to cut out the pads along the dotted lines. Outlines are actual size.
3. Place the paper cutouts on your face to determine if you prefer size Small or Large.



**Sleep Comfort Care Pad, Small**

Small pad should be 3.5" at its greatest width & 1.56" at its greatest height.

cut along dotted line



**Sleep Comfort Care Pad, Large**

Large pad should be 4.18" at its greatest width & 2.0" at its greatest height.

**Contact your CPAP mask provider to order today!**

**SEQUAL**<sup>®</sup>  
[www.sequal.com](http://www.sequal.com)